oligo

April 20, 2011

basecontent

Sequence Base Contents

Description

Function to compute the amounts of each nucleotide in a sequence.

Usage

basecontent (seq)

Arguments

seq

character vector of length n containg a valid sequence (A/T/C/G)

Value

matrix with n rows and 4 columns with the counts for each base.

Examples

```
sequences <- c("ATATATCCCCG", "TTTCCGAGC")
basecontent(sequences)</pre>
```

basicRMA

Simplified Interface to RMA

Description

Simple interface to RMA.

Usage

basicRMA(pmMat, pnVec, normalize = TRUE, background = TRUE, bgversion = 2, destr

2 boxplot-methods

Arguments

pmMat Matrix of intensities to be processed.

pnVec Probeset names.

normalize Logical flag: normalize?

background Logical flag: background adjustment?

bgversion Version of background correction.

destructive Logical flag: use destructive methods?

verbose Logical flag: verbose.
... Not currently used.

Value

Matrix.

Examples

```
set.seed(1)
pms <- matrix(rnorm(1000), nc=20)
colnames(pms) <- paste("sample", 1:20, sep="")
pns <- rep(letters[1:10], each=5)
res <- basicRMA(pms, pns, length(unique(pns)), TRUE, TRUE)
res[, 1:3]</pre>
```

boxplot-methods Boxplot

Description

Boxplot for observed (log-)intensities in a FeatureSet-like object (ExpressionFeatureSet, ExonFeatureSet, SnpFeatureSet, TilingFeatureSet) and ExpressionSet.

Usage

```
## S4 method for signature 'FeatureSet':
boxplot(x, which=c("pm", "mm", "bg", "both", "all"), transfo=log2, nsample=10000
## S4 method for signature 'ExpressionSet':
boxplot(x, which, transfo=identity, nsample=10000, ...)
```

Arguments

```
x a FeatureSet-like object or ExpressionSet object.
which character defining what probe types are to be used in the plot.
transfo a function to transform the data before plotting. See 'Details'.
nsample number of units to sample and build the plot.
... arguments to be passed to the default boxplot method.
```

chromosome 3

Details

The 'transfo' argument will set the transformation to be used. For raw data, 'transfo=log2' is a common practice. For summarized data (which are often in log2-scale), no transformation is needed (therefore 'transfo=identity').

Note

The boxplot methods for FeatureSet and Expression use a sample (via sample) of the probes/probesets to produce the plot. Therefore, the user interested in reproducibility is advised to use set.seed.

See Also

```
hist, image, sample, set.seed
```

chromosome

Accessor for chromosome information

Description

Returns chromosome information.

Usage

```
pmChr(object)
```

Arguments

object

TilingFeatureSet or SnpCallSet object

Details

chromosome () returns the chromosomal information for all probes and pmChr () subsets the output to the PM probes only (if a TilingFeatureSet object).

Value

Vector with chromosome information.

4 getX

darkColors

Create set of colors, interpolating through a set of preferred colors.

Description

Create set of colors, interpolating through a set of preferred colors.

Usage

```
darkColors(n)
seqColors(n)
```

Arguments

n

integer determining number of colors to be generated

Details

darkColors is based on the Dark2 palette in RColorBrewer, therefore useful to describe qualitative features of the data.

seqColors is based on Blues and generates a gradient of blues, therefore useful to describe quantitative features of the data.

Examples

```
x <- 1:10
y <- 1:10
cols1 <- darkColors(10)
cols2 <- seqColors(10)
plot(x, y, col=cols1, xlim=c(1, 11))
points(x+1, y, col=cols2)</pre>
```

getX

Accessors for physical array coordinates.

Description

Accessors for physical array coordinates.

Usage

```
getX(object, type)
getY(object, type)
```

Arguments

object FeatureSet object

type 'character' defining the type of the probes to be queried. Valid options are 'pm',

'mm', 'bg'

crlmm 5

Value

A vector with the requested coordinates.

Examples

```
## Not run:
x <- read.celfiles(list.celfiles())
theXpm <- getX(x, "pm")
theYpm <- getY(x, "pm")
## End(Not run)</pre>
```

crlmm

Genotype Calls

Description

Performs genotype calls via CRLMM (Corrected Robust Linear Model with Maximum-likelihood based distances).

Usage

```
crlmm(filenames, outdir, batch_size=40000, balance=1.5,
    minLLRforCalls=c(5, 1, 5), recalibrate=TRUE,
    verbose=TRUE, pkgname, reference=TRUE)
justCRLMM(filenames, batch_size = 40000, minLLRforCalls = c(5, 1, 5),
recalibrate = TRUE, balance = 1.5, phenoData = NULL, verbose = TRUE,
pkgname = NULL, tmpdir=tempdir())
```

Arguments

filenames character vector with the filenames.

outdir directory where the output (and some tmp files) files will be saved.
batch_size integer defining how many SNPs should be processed at a time.

recalibrate Logical - should recalibration be performed?

balance Control parameter to balance homozygotes and heterozygotes calls.

minLLRforCalls

Minimum thresholds for genotype calls.

verbose Logical.

phenoData phenoData object or NULL pkgname alt. pdInfo package to be used reference logical, defaulting to TRUE ...

tmpdir Directory where temporary files are going to be stored at.

Value

SnpCallSetPlus object.

6 getBaseProfile

```
getAffinitySplineCoefficients
```

Estimate affinity coefficients.

Description

Estimate affinity coefficients using sequence information and splines.

Usage

```
getAffinitySplineCoefficients(intensities, sequences)
```

Arguments

```
intensities Intensity matrix sequences Probe sequences
```

Value

Matrix with estimated coefficients.

See Also

getBaseProfile

```
getBaseProfile
```

Compute and plot nucleotide profile.

Description

Computes and, optionally, lots nucleotide profile, describing the sequence effect on intensities.

Usage

```
getBaseProfile(coefs, probeLength = 25, plot = FALSE, ...)
```

Arguments

```
coefs affinity spline coefficients.

probeLength length of probes
```

plot logical. Plots profile?

arguments to be passed to matplot.

Value

Invisibly returns a matrix with estimated effects.

getContainer 7

getContainer

Get container information for NimbleGen Tiling Arrays.

Description

Get container information for NimbleGen Tiling Arrays. This is useful for better identification of control probes.

Usage

```
getContainer(object, probeType)
```

Arguments

object A TilingFeatureSet or TilingFeatureSet object.

probeType String describing which probes to query ('pm', 'bg')

Value

'character' vector with container information.

getCrlmmSummaries Function to get CRLMM summaries saved to disk

Description

This will read the summaries written to disk and return them to the user as a SnpCallSetPlus or SnpCnvCallSetPlus object.

Usage

```
getCrlmmSummaries(tmpdir)
```

Arguments

tmpdir directory where CRLMM saved the results to.

Value

If the data were from SNP 5.0 or 6.0 arrays, the function will return a SnpCnvCallSetPlus object. It will return a SnpCallSetPlus object, otherwise.

getNgsColorsInfo

getNetAffx

NetAffx Biological Annotations

Description

Gets NetAffx Biological Annotations saved in the annotation package (Exon and Gene ST Affymetrix arrays).

Usage

```
getNetAffx(object, type = "probeset")
```

Arguments

object 'ExpressionSet' object (eg., result of rma())

type Either 'probeset' or 'transcript', depending on what type of summaries were

obtained.

Details

This retrieves NetAffx annotation saved in the (pd) annotation package - annotation(object). It is only available for Exon ST and Gene ST arrays.

The 'type' argument should match the summarization target used to generate 'object'. The 'rma' method allows for two targets: 'probeset' (target='probeset') and 'transcript' (target='core', target='full', target='extended').

Value

'AnnotatedDataFrame' that can be used as featureData(object)

Author(s)

Benilton Carvalho

getNgsColorsInfo

Helper function to extract color information for filenames on Nimble-Gen arrays.

Description

This function will (try to) extract the color information for NimbleGen arrays. This is useful when using read.xysfiles2 to parse XYS files for Tiling applications.

Usage

```
getNgsColorsInfo(path = ".", pattern1 = "_532", pattern2 = "_635", ...)
```

getPlatformDesign 9

Arguments

path path where to look for files

pattern1 pattern to match files supposed to go to the first channel

pattern2 pattern to match files supposed to go to the second channel

... extra arguments for list.xysfiles

Details

Many NimbleGen samples are identified following the pattern sampleID_532.XYS / sampleID_635.XYS. The function suggests sample names if all the filenames follow the standard above.

Value

A data.frame with, at least, two columns: 'channel1' and 'channel2'. A third column, 'sample-Names', is returned if the filenames follow the sampleID_532.XYS / sampleID_635.XYS standard.

Author(s)

Benilton Carvalho

bcarvalh@jhsph.edu>

```
getPlatformDesign Retrieve Platform Design object
```

Description

Retrieve platform design object.

Usage

```
getPlatformDesign(object)
getPD(object)
```

Arguments

object FeatureSet object

Details

Retrieve platform design object.

Value

```
platformDesign or PDInfo object.
```

10 image-methods

hist-methods	Density estimate
--------------	------------------

Description

Plot the density estimates for each sample

Usage

Arguments

```
x FeatureSet or ExpressionSet object
transfo a function to transform the data before plotting. See 'Details'.

nsample number of units to sample and build the plot.

which set of probes to be plotted ("pm", "mm", "bg", "both", "all").

... arguments to be passed to matplot
```

Details

The 'transfo' argument will set the transformation to be used. For raw data, 'transfo=log2' is a common practice. For summarized data (which are often in log2-scale), no transformation is needed (therefore 'transfo=identity').

Note

The hist methods for FeatureSet and Expression use a sample (via sample) of the probes/probesets to produce the plot (unless nsample > nrow(x)). Therefore, the user interested in reproducibility is advised to use set.seed.

image-methods	Display a pseudo-image of a microarray chip	
---------------	---	--

Description

Produces a pseudo-image (graphics::image) for each sample.

mmindex 11

Usage

```
## S4 method for signature 'FeatureSet':
image(x, which, transfo=log2, ...)
```

Arguments

X	FeatureSet object
which	integer indices of samples to be plotted (optional).
transfo	function to be applied to the data prior to plotting.
	parameters to be passed to image

mmindex

Accessors for PM, MM or background probes indices.

Description

Extracts the indexes for PM, MM or background probes.

Usage

```
mmindex(object, ...)
pmindex(object, ...)
bgindex(object, ...)
```

Arguments

```
object FeatureSet or DBPDInfo object
... Extra arguments, not yet implemented
```

Details

The indices are ordered by 'fid', i.e. they follow the order that the probes appear in the CEL/XYS files.

Value

A vector of integers representing the rows of the intensity matrix that correspond to PM, MM or background probes.

Examples

```
## How pm() works
## Not run:
x <- read.celfiles(list.celfiles())
pms0 <- pm(x)
pmi <- pmindex(x)
pms1 <- exprs(x)[pmi,]
identical(pms0, pms1)
## End(Not run)</pre>
```

12 justSNPRMA

mm

Accessors and replacement methods for the PM/MM/BG matrices.

Description

Accessors and replacement methods for the PM/MM/BG matrices.

Usage

```
mm(object, subset = NULL)
pm(object, subset = NULL, ...)
bg(object, subset = NULL)
mm(object) <-value
pm(object) <-value
bg(object) <-value</pre>
```

Arguments

```
object FeatureSet object.
subset Not implemented yet.
value matrix object.
... Extra arguments for future implementation.
```

Details

For all objects but TilingFeatureSet, these methods will return matrices. In case of TilingFeatureSet objects, the value is a 3-dimensional array (probes x samples x channels).

Examples

```
if (require(maqcExpression4plex) & require(pd.hg18.60mer.expr)) {
  xysPath <- system.file("extdata", package="maqcExpression4plex")
  xysFiles <- list.xysfiles(xysPath, full.name=TRUE)
  ngsExpressionFeatureSet <- read.xysfiles(xysFiles)
  pm(ngsExpressionFeatureSet)[1:10,]
}</pre>
```

justSNPRMA

Summarization of SNP data

Description

This function implements the SNPRMA method for summarization of SNP data. It works directly with the CEL files, saving memory.

Usage

```
justSNPRMA(filenames, verbose = TRUE, phenoData = NULL, normalizeToHapmap = TRUE
```

list.xysfiles 13

Arguments

filenames character vector with the filenames.

verbose logical flag for verbosity.

phenoData a phenoData object or NULL

normalizeToHapmap

Normalize to Hapmap? Should always be TRUE, but it's kept here for future use.

Value

SnpQSet or a SnpCnvQSet, depending on the array type.

Examples

```
## snprmaResults <- justSNPRMA(list.celfiles())</pre>
```

list.xysfiles

List XYS files

Description

Lists the XYS files.

Usage

```
list.xysfiles(...)
```

Arguments

... parameters to be passed to list.files

Details

The functions interface <code>list.files</code> and the user is asked to check that function for further details.

Value

Character vector with the filenames.

See Also

```
list.files
```

Examples

```
list.xysfiles()
```

14 oligo-package

Description

Create MA plots using a reference array (if one channel) or using channel2 as reference (if two channel).

Methods

object = "FeatureSet" ExpressionFeatureSet

oligo-package	The oligo package: a tool for low-level analysis of oligonucleotide arrays

Description

The **oligo** package provides tools to preprocess different oligonucleotide arrays types: expression, tiling, SNP and exon chips. The supported manufacturers are Affymetrix and NimbleGen.

It offers support to large datasets (when the **bigmemory** is loaded) and can execute preprocessing tasks in parallel (if, in addition to **bigmemory**, the **snow** package is also loaded).

Details

The package will read the raw intensity files (CEL for Affymetrix; XYS for NimbleGen) and allow the user to perform analyses starting at the feature-level.

Reading in the intensity files require the existence of data packages that contain the chip specific information (X/Y coordinates; feature types; sequence). These data packages packages are built using the **pdInfoBuilder** package.

For Affymetrix SNP arrays, users are asked to download the already built annotation packages from BioConductor. This is because these packages contain metadata that are not automatically created. The following annotation packages are available:

50 KXba - pd.mapping50 kxba.24050K Hind - pd.mapping50 khind.240250K Sty - pd.mapping250 k.sty250K Nsp - pd.mapping250 k.nsp GenomeWideSnp 5 (SNP 5.0) - pd.genomewidesnp.5 GenomeWideSnp 6 (SNP 6.0) - pd.genomewidesnp.6

For users interested in genotype calls for SNP 5.0 and 6.0 arrays, we strongly recommend the use use the **crlmm** package, which implements a more efficient version of CRLMM.

Author(s)

Benilton Carvalho - < carvalho@bclab.org>

References

Carvalho, B.; Bengtsson, H.; Speed, T. P. & Irizarry, R. A. Exploration, Normalization, and Genotype Calls of High Density Oligonucleotide SNP Array Data. Biostatistics, 2006.

plotM-methods 15

plotM-methods

Methods for Log-Ratio plotting

Description

The plotM methods are meant to plot log-ratios for different classes of data.

Methods

```
object = "SnpQSet", i = "character" Plot log-ratio for SNP data for sample i.
object = "SnpQSet", i = "integer" Plot log-ratio for SNP data for sample i.
object = "SnpQSet", i = "numeric" Plot log-ratio for SNP data for sample i.
object = "TilingQSet", i = "missing" Plot log-ratio for Tiling data for sample i.
```

pmAllele

Access the allele information for PM probes.

Description

Accessor to the allelic information for PM probes.

Usage

```
pmAllele(object)
```

Arguments

object

SnpFeatureSet or PDInfo object.

pmFragmentLength

Access the fragment length for PM probes.

Description

Accessor to the fragment length for PM probes.

Usage

```
pmFragmentLength(object)
```

Arguments

object

PDInfo object.

pmStrand

pmPosition

Accessor to position information

Description

pmPosition will return the genomic position for the (PM) probes.

Usage

```
pmPosition(object)
pmOffset(object)
```

Arguments

object

AffySNPPDInfo, TilingFeatureSet or SnpCallSet object

Details

```
pmPosition will return genomic position for PM probes on a tiling array. pmOffset will return the offset information for PM probes on SNP arrays.
```

pmStrand

Accessor to the strand information

Description

Returns the strand information on SNP arrays for PM probes (0 - sense / 1 - antisense).

Usage

```
pmStrand(object)
```

Arguments

object

AffySNPPDInfo object

summarize 17

summarize	Tools for microarray preprocessing	
-----------	------------------------------------	--

Description

Preprocess microarray data. Includes background correction, normalization and summarization methods.

Usage

```
backgroundCorrect(object, method="rma", copy=TRUE, verbose=TRUE, ...)
summarize(object, probes=rownames(object), method="medianpolish", verbose=TRUE,
normalize(object, method="quantile", copy=TRUE, verbose=TRUE, ...)
normalizeToTarget(object, target, method="quantile", copy=TRUE, verbose=TRUE)
```

Arguments

object	Object containing probe intensities to be preprocessed.
method	String determining which method to use at that preprocessing step.
target	Vector with the target distribution
probes	Character vector that identifies the name of the probes represented by the rows of object.
сору	Logical flag determining if data must be copied before processing (TRUE), or if data can be overwritten (FALSE).
verbose	Logical flag for verbosity.
	Arguments to be passed to methods.

Details

Number of rows of object must match the length of probes. Currently, only the following methods are implemented: - backgroundCorrectiong: rma.background - normalize: quantile - summarization: median-polish

Value

Expression matrix with length (unique (probes)) rows and ncol (object) columns.

Examples

```
ns <- 100
nps <- 1000
np <- 10
intensities <- matrix(rnorm(ns*nps*np, 8000, 400), nc=ns)
ids <- rep(as.character(1:nps), each=np)
bgCorrected <- backgroundCorrect(intensities)
normalized <- normalize(bgCorrected)
expression <- summarize(normalized, probes=ids)
intensities[1:20, 1:3]
expression[1:20, 1:3]
target <- rnorm(np*nps)
normalizedToTarget <- normalizeToTarget(intensities, target)</pre>
```

18 read.celfiles

probeNames Accessor to feature names

Description

Accessor to PM feature names.

Usage

```
probeNames(object, subset = NULL)
```

Arguments

object FeatureSet or DBPDInfo

subset not implemented yet.

read.celfiles Parser to CEL files

Description

Reads CEL files.

Usage

```
read.celfiles(..., filenames, pkgname, phenoData, featureData, experimentData, protocolData, notes, verbose=TRUE, sampleNames, rm.mask=FALSE, rm.outliers=FALSE, rm.extra=FALSE, checkType=TRUE)
```

read.celfiles2(channel1, channel2, pkgname, phenoData, featureData,
experimentData, protocolData, notes, verbose=TRUE, sampleNames,
rm.mask=FALSE, rm.outliers=FALSE, rm.extra=FALSE, checkType=TRUE)

Arguments

... names of files to be read.

filenames a character vector with the CEL filenames.

channel1 a character vector with the CEL filenames for the first 'channel' on a Tiling

application

channel2 a character vector with the CEL filenames for the second 'channel' on a

Tiling application

pkgname alternative data package to be loaded.

phenoData phenoData featureData

experimentData

experimentData

protocolData protocolData

readSummaries 19

```
notes

verbose logical

sampleNames character vector with sample names (usually better descriptors than the filenames)

rm.mask logical. Read masked?

rm.outliers logical. Remove outliers?

rm.extra logical. Remove extra?

checkType logical. Check type of each file? This can be time consuming.
```

Details

When using 'affyio' to read in CEL files, the user can read compressed CEL files (CEL.gz). Additionally, 'affyio' is much faster than 'affxparser'.

The function guesses which annotation package to use from the header of the CEL file. The user can also provide the name of the annotation package to be used (via the pkgname argument). If the annotation package cannot be loaded, the function returns an error. If the annotation package is not available from BioConductor, one can use the pdInfoBuilder package to build one.

Value

See Also

```
list.celfiles, read.xysfiles
```

Examples

```
if(require(pd.mapping50k.xba240) & require(hapmap100kxba)){
celPath <- system.file("celFiles", package="hapmap100kxba")
celFiles <- list.celfiles(celPath, full.name=TRUE)
affySnpFeatureSet <- read.celfiles(celFiles)
}</pre>
```

readSummaries

Read summaries generated by crlmm

Description

This function read the different summaries generated by crlmm.

20 read.xysfiles

Usage

```
readSummaries(type, tmpdir)
```

Arguments

type of summary of character class: 'alleleA', 'alleleB', 'alleleA-sense',

'alleleA-antisense', 'alleleB-sense', 'alleleB-antisense', 'calls', 'llr', 'conf'.

tmpdir directory containing the output saved by crlmm

Details

On the 50K and 250K arrays, given a SNP, there are probes on both strands (sense and antisense). For this reason, the options 'alleleA-sense', 'alleleA-antisense', 'alleleB-sense' and 'alleleB-antisense' should be used **only** with such arrays (XBA, HIND, NSP or STY).

On the SNP 5.0 and SNP 6.0 platforms, this distinction does not exist in terms of algorithm (note that the actual strand could be queried from the annotation package). For these arrays, options 'alleleA', 'alleleB' are the ones to be used.

The options calls, llr and conf will return, respectively, the CRLMM calls, log-likelihood ratios (for devel purpose **only**) and CRLMM confidence calls matrices.

Value

Matrix with values of summaries.

read.xysfiles Parser to XYS files

Description

NimbleGen provides XYS files which are read by this function.

Usage

```
read.xysfiles(..., filenames, pkgname, phenoData, featureData,
experimentData, protocolData, notes, verbose=TRUE, sampleNames,
checkType=TRUE)
```

read.xysfiles2(channel1, channel2, pkgname, phenoData, featureData,
experimentData, protocolData, notes, verbose=TRUE, sampleNames,
checkType=TRUE)

Arguments

	file names
filenames	character vector with filenames.
channel1	a character vector with the XYS filenames for the first 'channel' on a Tiling application $% \left(1\right) =\left(1\right) \left($
channel2	a character vector with the XYS filenames for the second 'channel' on a Tiling application

rma-methods 21

```
pkqname
                character vector with alternative PD Info package name
phenoData
               phenoData
featureData
               featureData
experimentData
               experimentData
protocolData protocolData
notes
               notes
verbose
               verbose
sampleNames
               character vector with sample names (usually better descriptors than the file-
               names)
               logical. Check type of each file? This can be time consuming.
checkType
```

Details

The function will read the XYS files provided by NimbleGen Systems and return an object of class FeatureSet.

The function guesses which annotation package to use from the header of the XYS file. The user can also provide the name of the annotation package to be used (via the pkgname argument). If the annotation package cannot be loaded, the function returns an error. If the annotation package is not available from BioConductor, one can use the pdInfoBuilder package to build one.

Value

See Also

```
list.xysfiles, read.celfiles
```

Examples

```
if (require(maqcExpression4plex) & require(pd.hg18.60mer.expr)) {
    xysPath <- system.file("extdata", package="maqcExpression4plex")
    xysFiles <- list.xysfiles(xysPath, full.name=TRUE)
    ngsExpressionFeatureSet <- read.xysfiles(xysFiles)
}</pre>
```

rma-methods

RMA - Robust Multichip Average algorithm

Description

Robust Multichip Average preprocessing methodology. This strategy allows background subtraction, quantile normalization and summarization (via median-polish).

22 rma-methods

Usage

```
## S4 method for signature 'ExonFeatureSet':
rma(object, background=TRUE, normalize=TRUE, subset=NULL, target="core")
## S4 method for signature 'ExpressionFeatureSet':
rma(object, background=TRUE, normalize=TRUE, subset=NULL)
## S4 method for signature 'GeneFeatureSet':
rma(object, background=TRUE, normalize=TRUE, subset=NULL, target="core")
## S4 method for signature 'SnpCnvFeatureSet':
rma(object, background=TRUE, normalize=TRUE, subset=NULL)
```

Arguments

object Exon/Expression/Gene/SnpCnv-FeatureSet object.
background Logical - perform RMA background correction?

normalize Logical - perform quantile normalization?

subset To be implemented.

target Level of summarization (only for Exon/Gene arrays)

Methods

signature (object = "ExonFeatureSet") When applied to an ExonFeatureSet object, rma can produce summaries at different levels: probeset (as defined in the PGF), core genes (as defined in the core.mps file), full genes (as defined in the full.mps file) or extended genes (as defined in the extended.mps file). To determine the level for summarization, use the target argument.

signature (object = "ExpressionFeatureSet") When used on an ExpressionFeatureSet object, rma produces summaries at the probeset level (as defined in the CDF or NDF files, depending on the manufacturer).

signature(object = "GeneFeatureSet") When applied to a GeneFeatureSet object, rma can produce summaries at different levels: probeset (as defined in the PGF) and
'core genes' (as defined in the core.mps file). To determine the level for summarization, use
the target argument.

signature (object = "SnpCnvFeatureSet") If used on a SnpCnvFeatureSet object (ie., SNP 5.0 or SNP 6.0 arrays), rma will produce summaries for the CNV probes. Note that this is an experimental feature for internal (and quick) assessment of CNV probes. We recommend the use of the 'crlmm' package, which contains a Copy Number tool specifically designed for these data.

References

Rafael. A. Irizarry, Benjamin M. Bolstad, Francois Collin, Leslie M. Cope, Bridget Hobbs and Terence P. Speed (2003), Summaries of Affymetrix GeneChip probe level data Nucleic Acids Research 31(4):e15

Bolstad, B.M., Irizarry R. A., Astrand M., and Speed, T.P. (2003), A Comparison of Normalization Methods for High Density O ligonucleotide Array Data Based on Bias and Variance. Bioinformatics 19(2):185-193

Irizarry, RA, Hobbs, B, Collin, F, Beazer-Barclay, YD, Antonellis, KJ, Scherf, U, Speed, TP (2003) Exploration, Normalizati on, and Summaries of High Density Oligonucleotide Array Probe Level Data. Biostatistics. Vol. 4, Number 2: 249-264

runDate 23

See Also

```
snprma
```

Examples

```
if (require(maqcExpression4plex) & require(pd.hg18.60mer.expr)) {
    xysPath <- system.file("extdata", package="maqcExpression4plex")
    xysFiles <- list.xysfiles(xysPath, full.name=TRUE)
    ngsExpressionFeatureSet <- read.xysfiles(xysFiles)
    summarized <- rma(ngsExpressionFeatureSet)
    show(summarized)
}</pre>
```

runDate

Date of scan

Description

Retrieves date information in CEL/XYS files.

Usage

```
runDate(object)
```

Arguments

object

'FeatureSet' object.

```
sequenceDesignMatrix
```

Create design matrix for sequences

Description

Creates design matrix for sequences.

Usage

```
sequenceDesignMatrix(seqs)
```

Arguments

seqs

character vector of 25-mers.

Details

This assumes all sequences are 25bp long.

The design matrix is often used when the objective is to adjust intensities by sequence.

24 snprma

Value

Matrix with length(seqs) rows and 75 columns.

Examples

```
genSequence <- function(x)
    paste(sample(c("A", "T", "C", "G"), 25, rep=TRUE), collapse="", sep="")
seqs <- sapply(1:10, genSequence)
X <- sequenceDesignMatrix(seqs)
Y <- rnorm(10, mean=12, sd=2)
Ydemean <- Y-mean(Y)
X[1:10, 1:3]
fit <- lm(Ydemean~X)
coef(fit)</pre>
```

mmSequence

Probe Sequeces

Description

Accessor to the (PM/MM/background) probe sequences.

Usage

```
mmSequence(object)
pmSequence(object, ...)
bgSequence(object, ...)
```

Arguments

```
object FeatureSet, AffySNPPDInfo or DBPDInfo object ... additional arguments
```

Value

A DNAStringSet containing the PM/MM/background probe sequence associated to the array.

snprma

Preprocessing SNP Arrays

Description

This function preprocess SNP arrays.

Usage

```
snprma(object, verbose = TRUE, normalizeToHapmap = TRUE)
```

snprma 25

Arguments

object SnpFeatureSet object verbose Verbosity flag. logical

normalizeToHapmap internal

Value

A SnpQSet object.

Index

*Topic IO	plotM-methods, 15
read.celfiles, 18	rma-methods, 21
read.xysfiles, 20	*Topic package
*Topic classif	oligo-package, 14
crlmm, 5	
getNetAffx, 8	backgroundCorrect (summarize), 17
runDate, 23	backgroundCorrect, ExpressionFeatureSet-method
*Topic file	(summarize), 17
list.xysfiles, 13	backgroundCorrect, FeatureSet-method
*Topic hplot	(summarize), 17
boxplot-methods, 2	backgroundCorrect, ff_matrix-method
darkColors,4	(summarize), 17
hist-methods, 10	backgroundCorrect, matrix-method
image-methods, 10	(summarize), 17
*Topic manip	backgroundCorrect-methods
basecontent, 1	(summarize), 17
basicRMA, 1	basecontent, 1
chromosome, 3	basicRMA, 1
getAffinitySplineCoefficients,	bg (mm), 12
6	bg, FeatureSet-method (mm), 12
<pre>getBaseProfile,6</pre>	bg, TilingFeatureSet-method (mm),
getContainer,7	12
<pre>getCrlmmSummaries,7</pre>	bg<- (mm), 12
${ t getNgsColorsInfo, 8}$	bg<-,FeatureSet,matrix-method
${ t getPlatformDesign}, { t 9}$	(mm), 12
getX,4	bg<-,TilingFeatureSet,array-method
justSNPRMA, 12	(mm), 12
mm, 12	bgindex(mmindex), 11
mmindex, 11	bgindex,DBPDInfo-method
mmSequence, 24	(mmindex), 11
pmAllele, 15	bgindex, FeatureSet-method
pmFragmentLength, 15	(mmindex), 11
pmPosition, 16	bgSequence (mmSequence), 24
pmStrand, 16	bgSequence,DBPDInfo-method
probeNames, 18	(mmSequence), 24
readSummaries, 19	bgSequence, ExonFeatureSet-method
sequenceDesignMatrix,23	(mmSequence), 24
snprma, 24	bgSequence, FeatureSet-method
summarize, 17	(mmSequence), 24
*Topic methods	bgSequence,GeneFeatureSet-method
boxplot-methods, 2	(mmSequence), 24
hist-methods, 10	boxplot,ExpressionSet-method
MAplot-methods, 14	(boxplot-methods), 2

INDEX 27

boxplot,FeatureSet-method	list.celfiles, 19
(boxplot-methods), 2	list.files, 13
boxplot-methods, 2	list.xysfiles, 13, 21
chromosome, 3	MAplot (MAplot-methods), 14
chromosome <- (chromosome), 3	MAplot, FeatureSet-method
chromosome <, Annotated Data Frame, chara	acter-methandlot-methods), 14
(chromosome), 3	MAplot-methods, 14
cleanPlatformName	mm, 12
(read.celfiles), 18	mm, FeatureSet-method (mm), 12
crlmm, 5	mm, TilingFeatureSet-method (mm),
darkColors,4	mm < - (mm), 12
1766' '1 0 1' 0 66' '1 1 6	mm<-, FeatureSet, matrix-method
getAffinitySplineCoefficients,6	(mm), 12
getBaseProfile,6	mm<-, TilingFeatureSet, array-method
getContainer,7	(mm), 12
getContainer, TilingFeatureSet-method	mmindex, 11
(getContainer), 7	mmindex, DBPDInfo-method
getContainer-methods	(mmindex), 11
(getContainer),7	mmindex, FeatureSet-method
getCrlmmSummaries,7	(mmindex), 11
getNetAffx, 8	mmSequence, 24
<pre>getNetAffx,ExpressionSet-method</pre>	mmSequence, AffySNPPDInfo-method
(getNetAffx), 8	
<pre>getNetAffx-methods (getNetAffx), 8</pre>	(mmSequence), 24
getNgsColorsInfo,8	mmSequence, DBPDInfo-method
getPD(getPlatformDesign), 9	(mmSequence), 24
<pre>getPlatformDesign, 9</pre>	mmSequence, FeatureSet-method
getPlatformDesign,FeatureSet-method	(mmSequence), 24
(getPlatformDesign), 9	7.1 () 17
getX,4	normalize (summarize), 17
getX, DBPDInfo-method (getX), 4	normalize, ff_matrix-method
getX, FeatureSet-method (getX), 4	(summarize), 17
getX-methods (getX), 4	normalize, matrix-method
getY (getX), 4	(summarize), 17
getY, DBPDInfo-method (getX), 4	normalize-methods (summarize), 17
getY, FeatureSet-method (getX), 4	normalizeToTarget(summarize),17
getY-methods (getX), 4	normalizeToTarget,ff_matrix-method
,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	(summarize), 17
hist, 3	normalizeToTarget, matrix-method
hist, ExpressionSet-method	(summarize),17
(hist-methods), 10	normalizeToTarget-methods
hist, FeatureSet-method	(summarize), 17
(hist-methods), 10	
hist-methods, 10	oligo-package, 14
image, 3	plotM(plotM-methods), 15
image, FeatureSet-method	plotM, SnpQSet, character-method
(image-methods), 10	(plotM-methods), 15
image-methods, 10	plotM, SnpQSet, integer-method
Image meerious, 10	(plotM-methods), 15
justCRLMM(crlmm), 5	plotM, SnpQSet, numeric-method
justSNPRMA, 12	(plotM-methods), 15

28 INDEX

plotM, TilingQSet, missing-method	pmStrand, AffySNPPDInfo-method
(plotM-methods), 15	(pmStrand), 16
plotM-methods, 15	probeNames, 18
pm (<i>mm</i>), 12	probeNames, DBPDInfo-method
pm, FeatureSet-method (mm), 12	(probeNames), 18
<pre>pm, TilingFeatureSet-method(mm),</pre>	<pre>probeNames, ExonFeatureSet-method</pre>
pm < -(mm), 12	<pre>probeNames,FeatureSet-method</pre>
pm<-,FeatureSet,matrix-method	(probeNames), 18
(mm), 12	<pre>probeNames,GeneFeatureSet-method</pre>
pm<-,TilingFeatureSet,array-method	(probeNames), 18
(mm), 12	
pmAllele, 15	read.celfiles, 18, 21
pmAllele, AffySNPPDInfo-method	read.celfiles2(read.celfiles), 18
(pmAllele), 15	read.xysfiles, 19, 20
·-	read.xysfiles2(read.xysfiles), 20
pmAllele, SnpFeatureSet-method	readSummaries, 19
(pmAllele), 15	rma (rma-methods), 21
pmChr (chromosome), 3	rma, ExonFeatureSet-method
pmChr, ExonFeatureSet-method	(rma-methods), 21
(chromosome), 3	rma, ExpressionFeatureSet-method
pmChr, FeatureSet-method	(rma-methods), 21
(chromosome), 3	rma, GeneFeatureSet-method
pmChr, GeneFeatureSet-method	(rma-methods), 21
(chromosome), 3	rma, SnpCnvFeatureSet-method
pmFragmentLength, 15	(rma-methods), 21
${\tt pmFragmentLength, AffySNPPDInfo-method}$	rma-methods, 21
(pmFragmentLength), 15	runDate, 23
pmindex (mmindex), 11	runDate, FeatureSet-method
pmindex, DBPDInfo-method	(runDate), 23
(mmindex), 11	runDate-methods (runDate), 23
pmindex, FeatureSet-method	Tunbace meenous (Tunbace), 25
(mmindex), 11	sample, 3
pmOffset (pmPosition), 16	seqColors (darkColors), 4
pmOffset, AffySNPPDInfo-method	sequenceDesignMatrix, 23
(pmPosition), 16	set . seed , 3
pmPosition, 16	snprma, 23, 24
pmPosition, ExpressionPDInfo-method	summarize, 17
(pmPosition), 16	summarize, ff_matrix-method
pmPosition, FeatureSet-method	(summarize), 17
(pmPosition), 16	summarize, matrix-method
pmPosition, TilingFeatureSet-method	(summarize), 17
(pmPosition), 16	summarize-methods (summarize), 17
pmPosition, TilingPDInfo-method	Summarize meenous (Summarize), 17
(pmPosition), 16	
· · · · · · · · · · · · · · · · · · ·	
pmSequence (mmSequence), 24	
pmSequence, AffySNPPDInfo-method	
(mmSequence), 24	
pmSequence, DBPDInfo-method	
(mmSequence), 24	
pmSequence, FeatureSet-method	
(mmSequence), 24	
pmStrand, 16	