Package 'scperturbR'

March 1, 2023

Type Package

Title E-Statistics for Seurat Objects

Version 0.1.0

Author Stefan Peidli

Maintainer Stefan Peidli <stefanpeidli@gmail.com>

Description R version of 'scperturb' tool for single-cell perturbation analysis. Contains wrappers for performing E-statistics for Seurat objects. More details on the method can be found in Peidli et al. (2023) <doi:10.1101/2022.08.20.504663> and in Székely and Rizzo (2004).

License MIT + file LICENSE

URL https://github.com/sanderlab/scPerturb/tree/master/package_r

BugReports https://github.com/sanderlab/scPerturb/issues

Depends R (>= 3.5.0) **Imports** Seurat, dplyr, rdist, energy, stats **Encoding** UTF-8 **Language** en-US

RoxygenNote 7.2.3

Suggests testthat (>= 3.0.0)

Config/testthat/edition 3

NeedsCompilation no

Repository CRAN

Date/Publication 2023-03-01 20:10:02 UTC

R topics documented:

edist	 	
etest	 	3

4

Index

edist

Description

Computes pairwise E-distances on a Seurat object. Computes E-distance between all groups in a Seurat object in space given by reduction.

Usage

```
edist(
  seurat_object,
  groupby = "perturbation",
  reduction = "pca",
  sample_correction = FALSE,
  verbose = TRUE
)
```

Arguments

seurat_object	An object of class Seurat.	
groupby	An object of class character. Points to the column in the Seurat object's meta data that contains the group labels.	
reduction	An object of class character. The reduction / embedding in seurat_object that is used to compute the E-distance in.	
sample_correction		
	An object of class logical. If TRUE, the E-distances are corrected for sample size. Will make it not a proper distance, leads to negative values.	
verbose	An object of class logical. If TRUE, prints messages. Default is TRUE.	

Value

Returns an object of class data.frame. For each group contains the E-test p-value and the E-distance to control group.

Examples

Add some code illustrating how to use the function

etest

Description

Performs E-testing on a Seurat object. Computes E-test statistics for each group in a Seurat object, using the E-distance in space given by reduction to the group defined by control.

Usage

```
etest(
  seurat_object,
  groupby = "perturbation",
  control = "control",
  reduction = "pca",
  verbose = TRUE,
  permutations = 1000
)
```

Arguments

seurat_object	An object of class Seurat.
groupby	An object of class character. Points to the column in the Seurat object's meta data that contains the group labels.
control	An object of class character. The group that is used as the control.
reduction	An object of class character. The reduction / embedding in seurat_object that is used to compute the E-distance in.
verbose	An object of class logical. If TRUE, prints messages. Default is TRUE.
permutations	An object of class integer. The number of permutations used to compute the p-value. Default is 1000.

Value

Returns an object of class data.frame. For each group contains the E-test p-value and the E-distance to control group.

Examples

Add some code illustrating how to use the function

Index

edist,2 etest,3