Package 'alabaster.matrix'

July 23, 2025

Title Load and Save Artifacts from File

Version 1.8.0

Date 2025-03-04

License MIT + file LICENSE

Description

Save matrices, arrays and similar objects into file artifacts, and load them back into memory.

This is a more portable alternative to serialization of such objects into RDS files.

Each artifact is associated with metadata for further interpretation;

downstream applications can enrich this metadata with context-specific properties.

Depends alabaster.base

Imports methods, BiocGenerics, S4Vectors, DelayedArray (>= 0.33.3),

S4Arrays, SparseArray (>= 1.5.22), rhdf5 (>= 2.47.1),

HDF5Array, Matrix, Rcpp

Suggests testthat, knitr, BiocStyle, chihaya, BiocSingular,

ResidualMatrix

LinkingTo Rcpp

VignetteBuilder knitr

RoxygenNote 7.3.2

Encoding UTF-8

biocViews DataImport, DataRepresentation

git_url https://git.bioconductor.org/packages/alabaster.matrix

git_branch RELEASE_3_21

git_last_commit 2c0db5d

git_last_commit_date 2025-04-15

Repository Bioconductor 3.21

Date/Publication 2025-07-23

Author Aaron Lun [aut, cre]

Maintainer Aaron Lun <infinite.monkeys.with.keyboards@gmail.com>

2 AmalgamatedArray

Contents

Amal	gamatedArray	Amalgamated array class	
Index			22
	writeSparseMatrix.		. 20
	_		
	saveDelayedArray .		. 13
	saveArray		. 11
	ReloadedArraySeed		. 10
	recycleHdf5Files .		. 9
	readSparseMatrix .		. 8
	•	erations	
	-		
		1	
	AmalgamatedArray		. 2

Description

Implements an amalgamated array, equivalent to a delayed combination of DelayedArray objects. It allows stageObject to save a combination of multiple matrices without actually aggregating their data into a single file.

Constructors

AmalgamatedArraySeed(..., along=1) accepts any number of named array-like objects and returns a AmalgamatedArraySeed. Each object corresponds to a block and should be named accordingly; names should be unique and non-empty. The along argument specifies the dimension in which matrices should be combined - for matrices, this is 1 for rows, 2 for columns.

AmalgamatedArray(..., along=1) accepts any number of named array-like objects and returns a AmalgamatedArray. Alternatively, a single AmalgamatedArraySeed may be provided in

Functions

componentNames(x) will return a character vector of names of component arrays in a AmalgamatedArray(Seed) object x.

extractComponents(x) will return a named list of array-like objects, corresponding to the component arrays used to construct the AmalgamatedArray(Seed) object x.

stageObject(x, dir, path, child = FALSE) will save the AmalgamatedArray x and its components into the path inside dir. Each component array is staged into its own subdirectory inside path. createRawArraySeed 3

Comments on usage

The AmalgamatedArraySeed is closely related to (and in fact, is a subclass of) the DelayedAbind class. This means that we can leverage many of the DelayedArray methods for handling the delayed bind. In theory, we could just use a DelayedAbind directly and save it with chihaya in stageObject (via preserveDelayedOperations(TRUE)). However, this provides fewer opportunities for tracking and manipulating the samples. It also saves the per-sample matrices into a single file, which eliminates possibilities for per-file deduplication and linking, e.g., with recycleHdf5Files(TRUE).

Author(s)

Aaron Lun

Examples

```
first <- Matrix::rsparsematrix(10, 10, 0.1)
second <- Matrix::rsparsematrix(10, 20, 0.1)
mat <- AmalgamatedArray(list(foo = first, bar = second), along=2)
mat

componentNames(mat)
out <- extractComponents(mat)
lapply(out, dim)</pre>
```

createRawArraySeed

Array loading utilities

Description

Utilities for loading an array saved by stageObject.

Usage

```
.createRawArraySeed(info, project, names = TRUE)
.extractArrayDimnames(path, group, ndim)
```

Arguments

info	A named list of metadata for this array.
project	Any argument accepted by the acquisition functions, see ?acquireFile. By default, this should be a string containing the path to a staging directory.
names	Logical scalar indicating whether the seed should be annotated with dimnames (if available).
path	String containing the path to the file containing said array.
group	String containing the name of the group with the dimnames.
ndim	Integer scalar specifying the number of dimensions.

4 DelayedMask

Details

For .createArraySeed, the array should be one of:

- hdf5_dense_array
- hdf5_sparse_matrix
- hdf5_delayed_array
- amalgamated_array

For delayed arrays, the file may contain a seed array with the "custom alabaster local array" type. This should have a path dataset containing a relative path to another array in the same project, which is loaded and used as the seed for this delayed array. Callers can overwrite this behavior by setting "custom alabaster local array" in the knownArrays from **chihaya** before calling .createRawArraySeed.

For .extractArrayDimnames, path is expected to be a HDF5 file with a group specified by group. Each child of this group is a string dataset named after a (0-indexed) dimension, containing the names for that dimension.

Value

- . createRawArraySeed returns a seed that can be used in the DelayedArray constructor.
- . extractArrayDimnames returns a list of character vectors or NULL, containing the dimnames.

Author(s)

Aaron Lun

Examples

```
# Staging an array as an example:
dir <- tempfile()
dir.create(dir)
mat <- array(rpois(10000, 10), c(50, 20, 10))
meta <- stageObject(mat, dir, "whee")
# Loading it back as a DelayedArray seed:
.createRawArraySeed(meta, project=dir)</pre>
```

DelayedMask

Delayed masking

Description

Delayed masking of missing values, based on replacement of placeholder values with NA. This allows missingness to be encoded in frameworks without the same concept of NA as R.

Usage

```
DelayedMask(x, placeholder)
```

Arguments

x An existing **DelayedArray** seed.

placeholder Placeholder value to replace with NA. This should be of the same type as type(x).

Details

If is.na(placeholder) is true for double-precision x, masking is performed for all values of x where is.na is true. This includes both NaNs and NAs; no attempt is made to distinguish between the NaN payloads.

Currently, an error is raised for any integer x that produces non-missing values of -2^31 without a placeholder of NA_integer_. This is because R cannot distinguish the integer -2^31 from an integer-type NA.

Value

A DelayedMask object, to be wrapped in a DelayedArray.

Author(s)

Aaron Lun

Examples

```
original <- DelayedArray(matrix(rpois(40, lambda=2), ncol=5))
original
masked <- DelayedMask(original, 0)
DelayedArray(masked)</pre>
```

preserveDelayedOperations

Preserve delayed operations during staging

Description

Preserve delayed operations via chihaya when staging a DelayedArray with stageObject.

Usage

```
preserveDelayedOperations(preserve)
```

Arguments

preserve

Whether to preserve delayed operations using the **chihaya** specification.

6 readArray

Details

By default, any DelayedArray in stageObject will be saved as a new dense array or sparse matrix. However, if this option is enabled, DelayedArrays will instead be saved in the **chihaya** specification, where the delayed operations are themselves stored in the HDF5 file (see https://artifactdb.github.io/chihaya/ for details).

The **chihaya** specification is more complicated to parse but can be helpful in reducing disk usage. One simple example is to avoid sparsity-breaking or integer-to-float operations by storing their delayed representations in the file. If the seed matrix is derived from some immutable reference location, advanced users can even store links to that location instead of duplicating the seed data.

Value

Logical scalar indicating whether delayed operations are to be preserved by the DelayedArray method. If preserve is supplied, it is used to set this scalar, and the *previous* value of the scalar is invisibly returned.

Author(s)

Aaron Lun

Examples

```
preserveDelayedOperations()
old <- preserveDelayedOperations(TRUE)
preserveDelayedOperations()
preserveDelayedOperations(old)</pre>
```

readArray

Read a dense array from disk

Description

Read a dense high-dimensional array from its on-disk representation. This is usually not directly called by users, but is instead called by dispatch in readObject.

Usage

```
readArray(path, metadata, ...)
```

Arguments

path String containing a path to a directory, itself created by the saveObject method

for a dense array.

metadata Named list of metadata for this object, see readObject for more details.

... Further arguments, ignored.

readDelayedArray 7

Value

A dense file-backed ReloadedArray.

Author(s)

Aaron Lun

See Also

"saveObject, array-method", to create the directory and its contents.

Examples

```
arr <- array(rpois(10000, 10), c(50, 20, 10))
dimnames(arr) <- list(
   paste0("GENE_", seq_len(nrow(arr))),
   letters[1:20],
   NULL
)
dir <- tempfile()
saveObject(arr, dir)
readObject(dir)</pre>
```

readDelayedArray

Read a delayed array from disk

Description

Read a delayed high-dimensional array from its on-disk representation. This is usually not directly called by users, but is instead called by dispatch in readObject.

Usage

```
readDelayedArray(path, metadata, delayed_array.reload.args = list(), ...)
```

Arguments

path String containing a path to a directory, itself created by the saveObject method for a delayed array.

metadata Named list of metadata for this object, see readObject for more details.

delayed_array.reload.args

Named list of arguments to be passed to reloadDelayedObject.

... Further arguments, ignored.

8 readSparseMatrix

Value

A multi-dimensional array-like object.

Author(s)

Aaron Lun

See Also

```
"saveObject, DelayedArray-method", to create the directory and its contents. reloadDelayedObject, for the methods to reload each delayed operation.
```

Examples

```
arr <- array(rpois(10000, 10), c(50, 20, 10))
dimnames(arr) <- list(
   paste0("GENE_", seq_len(nrow(arr))),
   letters[1:20],
   NULL
)
dir <- tempfile()
saveObject(arr, dir)
readObject(dir)</pre>
```

readSparseMatrix

Read a sparse matrix from disk

Description

Read a sparse matrix from its on-disk representation. This is usually not directly called by users, but is instead called by dispatch in readObject.

Usage

```
readSparseMatrix(path, metadata, ...)
```

Arguments

path String containing a path to a directory, itself created by the saveObject method for a spars matrix.

metadata Named list of metadata for this object, see readObject for more details.

... Further arguments, ignored.

Value

A sparse ReloadedMatrix object.

recycleHdf5Files 9

Author(s)

Aaron Lun

See Also

"saveObject, sparseMatrix-method", to create the directory and its contents.

Examples

```
mat <- Matrix::rsparsematrix(100, 200, density=0.2)
rownames(mat) <- paste0("GENE_", seq_len(nrow(mat)))
dir <- tempfile()
saveObject(mat, dir)
readObject(dir)</pre>
```

recycleHdf5Files

Recycle existing HDF5 files

Description

Re-use existing files in HDF5-backed arrays rather than reserializing them in stageObject.

Usage

```
recycleHdf5Files(recycle)
```

Arguments

recycle

Whether to recycle existing files for HDF5-backed DelayedArrays.

Details

If this options is enabled, stageObject will attempt to link/copy existing files for any HDF5-backed DelayedArray instances - most specifically, HDF5Array objects and H5SparseMatrix objects using the 10X format. This avoids re-serialization of the data for faster staging. It also allows advanced users to add their own customizations into the HDF5 file during staging, as long as they do not interfere with loadArray.

By default, this option is disabled as the properties of the existing file are not known in the general case. In particular, the file might contain other groups/datasets that are irrelevant, and use up extra disk space if copied; or confidential, and should not be stored in the staging directory. Users should only enable this option if they have full control over the generation and contents of the backing HDF5 files.

Also note that any dimnames on x will be ignored during recycling.

10 ReloadedArraySeed

Value

Logical scalar indicating whether HDF5 files are to be reused. If recycle is supplied, it is used to set this scalar, and the *previous* value of the scalar is invisibly returned.

Author(s)

Aaron Lun

Examples

```
recycleHdf5Files()
old <- recycleHdf5Files(TRUE)
recycleHdf5Files()
recycleHdf5Files(old)</pre>
```

ReloadedArraySeed

Reloaded alabaster array

Description

An array that was reloaded from disk by the readObject function. This allows methods to refer to the existing on-disk representation by inspecting the path. For example, saveObject can just copy/link to the files instead of repeating the saving process.

Usage

```
ReloadedArraySeed(path, seed = NULL, ...)
ReloadedArray(path, seed = NULL, ...)
```

Arguments

path	String containing a path to the directory with the on-disk array representation. Alternatively an existing ReloadedArraySeed, which is returned without modification.
seed	Contents of the loaded array, e.g., as an ordinary R array, a DelayedArray or a sparse matrix. If NULL, this is obtained by calling readObject.
	Further arguments to pass to readObject when seed=NULL.

Details

The ReloadedArraySeed is a DelayedUnaryIsoOp subclass that will just forward all operations to the underlying seed. Its main purpose is to track the path that was originally used to generate seed, which enables optimizations for methods that need to operate on the files.

One obvious optimization is the specialization of saveObject on ReloadedArray instances. Instead of loading the array data back into the R session and saving it again, the saveObject method can just

saveArray 11

link or copy the existing files. This behavior is controlled by the ReloadedArray.reuse.files= option in the saveObject method, which can be any of the choices for action= in cloneDirectory.

It may also be "none" to ignore existing files and just save the contents by calling "saveObject, DelayedArray-method".

Value

For the constructors, an instance of the ReloadedArraySeed or ReloadedArray.

Examples

```
arr <- array(rpois(10000, 10), c(50, 20, 10))
dir <- tempfile()
saveObject(arr, dir)
obj <- readArray(dir)
obj
DelayedArray::showtree(obj)</pre>
```

saveArray

Save a multi-dimensional array to disk

Description

Save a high-dimensional array to its on-disk representations.

Usage

```
## S4 method for signature 'array'
saveObject(
    X,
    path,
    array.dedup.session = NULL,
    array.dedup.action = "link",
    array.character.vls = FALSE,
    ...
)

## S4 method for signature 'denseMatrix'
saveObject(
    x,
    path,
    array.dedup.session = NULL,
    array.dedup.action = "link",
    ...
)
```

12 saveArray

Arguments

x An integer, numeric, logical or character array. Alternatively, any of the dense-

Matrix subclasses from the Matrix package.

path String containing the path to a directory in which to save x.

array.dedup.session

A session object created by createDedupSession, specifying which objects should be deduplicated if the same x is encountered multiple times.

array.dedup.action

String specifying how deduplication should occur, see options for the action= argument in cloneDirectory.

array.character.vls

Logical scalar indicating whether to save character arrays in the custom variable length string (VLS) array format. If NULL, this is determined based on a comparison of the expected storage against a fixed length array.

. . . Further arguments, currently ignored.

Details

When saveObject is called multiple times on the same x, this method can avoid re-saving the array if a user supplies a deduplication session in array.dedup.session=. Instead, the method will link or copy (depending on the choice of array.dedup.action=) the files produced by the first saveObject call to the new path=. This saves time and reduces disk usage, and is particularly useful when saving complex data structures like a SummarizedExperiment. For example, if we have a SummarizedExperiment that contains multiple copies of the same array, we can pass a array.dedup.session= to the saveObject call on the SummarizedExperiment. This will instruct the internal saveObject calls to only write the array to disk once and subsequently make copies or links for all duplicates of that array. The same approach can be applied to deduplicate seeds in a DelayedArray, see ?"storeDelayedObject" for details.

Value

x is saved to path and NULL is invisibly returned.

Author(s)

Aaron Lun

See Also

readArray, to read the directory contents back into the R session.

Examples

```
mat <- array(rpois(10000, 10), c(50, 20, 10))
dimnames(mat) <- list(
   paste0("GENE_", seq_len(nrow(mat))),
   letters[1:20],
   NULL</pre>
```

saveDelayedArray 13

```
dir <- tempfile()
saveObject(mat, dir)
list.files(dir)</pre>
```

saveDelayedArray

Save DelayedArrays to disk

Description

Save DelayedArray objects to their on-disk representation.

Usage

```
## S4 method for signature 'DelayedArray'
saveObject(
    X,
    path,
    DelayedArray.dispatch.pristine = TRUE,
    DelayedArray.preserve.ops = FALSE,
    DelayedArray.force.external = FALSE,
    DelayedArray.store.args = list(),
    array.dedup.session = NULL,
    array.dedup.action = NULL,
    ...
)
```

Arguments

x A DelayedArray object.

path String containing a path to a directory in which to save x.

DelayedArray.dispatch.pristine

Logical scalar indicating whether to call the saveObject methods of seeds of pristine arrays.

DelayedArray.preserve.ops

Logical scalar indicating whether delayed operations should be preserved ondisk.

DelayedArray.force.external

Logical scalar indicating to save the seeds of x as external arrays. This is passed directly to storeDelayedObject as the save.external.array= argument, see ?storeDelayedObject for details.

DelayedArray.store.args

More named arguments to pass to storeDelayedObject.

14 saveSparseMatrix

array.dedup.session, array.dedup.action

Arguments controlling deduplication of x, see ?"saveObject, array-method" for details. If x is not a duplicate of an existing object, these arguments will be passed to further methods as described for

Further arguments passed to storeDelayedObject as external.save.args, if the delayed operations are to be preserved; otherwise, they are passed to saveObject,array-method or saveObject,sparseMatrix-method.

Details

Supplying array.dedup.session= by itself is only guaranteed to deduplicate x itself and may not deduplicate its seeds. Users should combine this with DelayedArray.force.external=TRUE to force seeds to be saved via saveObject, which exposes the seeds to the deduplication machinery in their respective saveObject methods. Check out ?"storeDelayedObject" for more details.

Value

x is saved to path and NULL is invisibly returned.

Author(s)

Aaron Lun

See Also

storeDelayedObject, for the methods to save each delayed operation.

Examples

```
mat <- Matrix::rsparsematrix(100, 200, density=0.2)
rownames(mat) <- paste0("GENE_", seq_len(nrow(mat)))
dmat <- DelayedArray::DelayedArray(mat) * 1

dir <- tempfile()
saveObject(dmat, dir, delayed.preserve.ops=TRUE)
list.files(dir)</pre>
```

saveSparseMatrix

Save a sparse matrix to disk

Description

Save a sparse matrix to its on-disk representations.

saveSparseMatrix 15

Usage

```
## S4 method for signature 'sparseMatrix'
saveObject(
    x,
    path,
    array.dedup.session = NULL,
    array.dedup.action = "link",
    ...
)

## S4 method for signature 'SVT_SparseMatrix'
saveObject(
    x,
    path,
    array.dedup.session = NULL,
    array.dedup.action = "link",
    ...
)
```

Arguments

for details.

... Further arguments, currently ignored.

Value

x is saved to path and NULL is invisibly returned.

Author(s)

Aaron Lun

See Also

readSparseMatrix, to read the directory contents back into the R session.

Examples

```
mat <- Matrix::rsparsematrix(100, 200, density=0.2)
rownames(mat) <- paste0("GENE_", seq_len(nrow(mat)))

dir <- tempfile()
saveObject(mat, dir)
list.files(dir)</pre>
```

16 storeDelayedObject

storeDelayedObject

Store/reload a DelayedArray

Description

Store or reload the delayed operations or array-like seeds of a DelayedArray in an existing HDF5 file.

Usage

```
storeDelayedObject(x, handle, name, ...)
reloadDelayedObject(handle, name, version = package_version("1.1"), ...)
reloadDelayedObjectFunctionRegistry(type = c("operation", "array"))
registerReloadDelayedObjectFunction(
  type = c("operation", "array"),
  subtype,
 fun,
  existing = c("old", "new", "error")
)
## S4 method for signature 'array'
storeDelayedObject(
 Х,
 handle,
 name,
 version = package_version("1.1"),
 save.external.array = FALSE,
)
## S4 method for signature 'ANY'
storeDelayedObject(
  Х,
 handle,
 name,
 version = package_version("1.1"),
 external.save.args = list(),
 external.dedup.session = NULL,
 external.dedup.action = NULL,
)
```

storeDelayedObject 17

```
altStoreDelayedObjectFunction(store)
altStoreDelayedObject(...)
altReloadDelayedObjectFunction(reload)
altReloadDelayedObject(...)
```

Arguments

x Any of the delayed operation/array classes from **DelayedArray**.

handle An **rhdf5** handle of a HDF5 file to save into (for storeDelayedObject) or load

from (for reloadDelayedObject).

name String containing the name of the group in file to save into (for storeDelayedObject)

or load from (for reloadDelayedObject).

.. For storeDelayedObject and reloadDelayedObject, additional arguments to

be passed to specific methods.

For altStoreDelayedObject and altReloadDelayedObject, arguments to be

passed to the alternative functions.

version Package version of the **chihava** format to use when storing or reloading delayed

objects. When reloading, the version should be retrieved from the attributes of

the outermost group, typically by readDelayedArray.

type String specifying the type of delayed object, i.e., operation or array. This corre-

sponds to delayed_type type in the **chihaya** attributes.

subtype String specifying the subtype of the delayed object, This corresponds to delayed_array

or delayed_operation type (depending on type) in the **chihaya** attributes.

fun Function to reload a delayed object. This should accept the same arguments as

reloadDelayedObject and should return a delayed array (if type="array") or operation (otherwise). It may also be NULL to delete an existing entry in the

registry.

existing Logical scalar indicating the action to take if a function has already been regis-

tered for type and subtype - keep the old or new function, or throw an error.

save.external.array

Logical scalar indicating whether to save an array-like seed as an external seed,

even if a dedicated storeDelayedObject method is available.

external.save.args

Named list of further arguments to pass to altSaveObject when saving an ex-

ternal seed.

external.dedup.session

Deprecated, set external.save.args\$array.dedup.session instead.

external.dedup.action

Deprecated, set external.save.args\$array.dedup.action instead.

store Function (typically a generic) to store delayed objects to file. This should accept

the same arguments as storeDelayedObject.

reload Function to reload delayed objects from file. This should accept the same argu-

ments as reloadDelayedObject.

18 storeDelayedObject

Value

For storeDelayedObject and altStoreDelayedObject, the contents of x are saved to file, and NULL is invisibly returned.

For reloadDelayedObject and altReloadDelayedObject, a delayed operation or DelayedArray is returned.

For altStoreDelayedObjectFunction, the current store function is returned if store is missing. Otherwise, store is set as the current store function and the previous store function is returned.

For altReloadDelayedObjectFunction, the current reload function is returned if reload is missing. Otherwise, reload is set as the current reload function and the previous reload function is returned.

Customization

Developers can easily extend **alabaster.matrix** to new delayed objects by writing new methods for storeDelayedObject. Methods should save the contents of the delayed object to the HDF5 file in the **chihaya** format. Each new store method typically requires a corresponding reloading function to be registered via registerReloadDelayedObjectFunction, so that reloadDelayedObject knows how to reconstitute the object from file.

Application developers can customize the process of storing/reloading delayed objects by specifying alternative functions in altReloadDelayedObjectFunction and altStoreDelayedObjectFunction. For example, if we want to preserve all delayed objects except for DelayedSubset, we could replace storeDelayedObject with an altStoreDelayedObject that realizes any DelayedSubset instance into an ordinary matrix. This is analogous to the overrides for altReadObject and altSaveObject.

Extension developers (i.e., those who write new methods for storeDelayedObject or new functions for reloadDelayedObject) should generally use altStoreDelayedObject and altReloadDelayedObject in their method/funcion bodies. This ensures that any custom overrides specified by application developers are still respected in the extensions to **alabaster.matrix**.

External seeds

Whenever storeDelayedObject encounters a delayed operation or array-like seed for which it has no methods, the ANY method will save the delayed object as an "external seed". The array is saved via altSaveObject into a seeds directory next to the file associated with handle. A reference to this external location is then stored in the name group inside handle.

Users can force this behavior for all array-like seeds by specifying save.external.array=TRUE. This instructs storeDelayedObject to save everything as external seeds, including those arrays for which it has methods. Doing so can be beneficial to enable deduplication, e.g., when two delayed arrays perform different operations on the same underlying seed. By saving the seeds externally, file management systems can identify the redundancy to save storage space.

Advanced users can explicitly deduplicate external seeds by setting both save.external.array=TRUE and supplying array.dedup.session= in external.save.args=. The array.dedup.session object is filled up with unique seeds as storeDelayedObject is called on various DelayedArrays (see ?"saveObject,array-method" for details). Whenever a duplicate seed is encountered, it is not saved again, but is instead linked or copied from the file path associated with the identical external seed. For example, a new session can be created when saving a SummarizedExperiment to deduplicate seeds across its assays.

WrapperArraySeed 19

When external seeds are encountered by reloadDelayedObject, they are loaded as ReloadedArrays (or some variant thereof) by altReadObject. Users can forcibly realize the reloaded seed into memory by passing custom.takane.reload=TRUE in ... for the reloadDelayedObject call. This is occasionally helpful for providing a more faithful roundtrip from file back into memory.

Author(s)

Aaron Lun

See Also

saveObject, DelayedArray-method and readDelayedArray, where these methods are used. https://artifactdb.github.io/chihaya/, for the file format specification of delayed objects.

Examples

```
library(DelayedArray)
X <- DelayedArray(matrix(runif(100), ncol=20))
Y <- cbind(X, DelayedArray::ConstantArray(value=50, c(5, 10)))
library(rhdf5)
temp <- tempfile()
dir.create(temp)

fpath <- file.path(temp, "foo.h5")
fhandle <- H5Fcreate(fpath)
storeDelayedObject(Y@seed, fhandle, "YAY")
rhdf5::h5ls(fhandle)
H5Fclose(fhandle)

fhandle <- H5Fopen(fpath, "H5F_ACC_RDONLY")
reloadDelayedObject(fhandle, "YAY")
H5Fclose(fhandle)</pre>
```

WrapperArraySeed

DelayedArray wrapper seed

Description

Deprecated, use DelayedUnaryIsoOp instead.

20 writeSparseMatrix

writeSparseMatrix

Write a sparse matrix

Description

Writes a sparse matrix to file in a compressed sparse format.

Usage

```
writeSparseMatrix(
    X,
    file,
    name,
    chunk = 10000,
    column = TRUE,
    tenx = FALSE,
    guess.integer = TRUE
)
```

Arguments

X	A sparse matrix of some sort. This includes sparse DelayedMatrix objects.
file	String containing a path to the HDF5 file. The file is created if it is not already present.
name	String containing the name of the group to store x.
chunk	Integer scalar specifying the chunk size for the indices and values.
column	Logical scalar indicating whether to store as compressed sparse column format.
tenx	Logical scalar indicating whether to use the $10\mathrm{X}$ compressed sparse column format.
guess.integer	Logical scalar specifying whether to guess an appropriate integer type from x.

Details

This writes a sparse matrix to file in various formats:

- column=TRUE and tenx=FALSE uses H5AD's csr_matrix format.
- column=FALSE and tenx=FALSE uses H5AD's csc_matrix format.
- tenx=TRUE uses 10X Genomics' HDF5 matrix format.

For the first two formats, the apparent transposition is deliberate, because columns in R are interpreted as rows in H5AD. This allows us to retain consistency the interpretation of samples (columns in R, rows in H5AD) and features (vice versa). Constructors for classes like H5SparseMatrix will automatically transpose so no extra work is required.

If guess.integer=TRUE, we attempt to save x's values into the smallest type that will accommodate all of its values. If x only contains unsigned integers, we will attempt to save either 8-, 16- or 32-bit

writeSparseMatrix 21

unsigned integers. If x contains signed integers, we will fall back to 32-bit signed integers. For all other values, we will fall back to double-precision floating point values.

We attempt to save x's indices to unsigned 16-bit integers if the relevant dimension of x is small enough. Otherwise we will save it as an unsigned 32-bit integer.

Value

A NULL invisibly. The contents of x are written to name in file.

Author(s)

Aaron Lun

Examples

```
library(Matrix)
x <- rsparsematrix(100, 20, 0.5)
tmp <- tempfile(fileext=".h5")
writeSparseMatrix(x, tmp, "csc_matrix")
writeSparseMatrix(x, tmp, "csr_matrix", column=FALSE)
writeSparseMatrix(x, tmp, "tenx_matrix", tenx = TRUE)

rhdf5::h5ls(tmp)
library(HDF5Array)
H5SparseMatrix(tmp, "csc_matrix")
H5SparseMatrix(tmp, "csr_matrix")
H5SparseMatrix(tmp, "tenx_matrix")</pre>
```

Index

.createRawArraySeed	createExternalSeedDedupSession	
(createRawArraySeed), 3	(storeDelayedObject), 16	
.extractArrayDimnames	createRawArraySeed, 3	
(createRawArraySeed), 3		
	DelayedAbind, 3	
acquireFile, 3	DelayedArray, 4, 5, 10, 13, 16, 18	
altReadObject, 18, 19	DelayedArray,AmalgamatedArraySeed-method	
altReloadDelayedObject	(AmalgamatedArray), 2	
(storeDelayedObject), 16	DelayedArray,ReloadedArraySeed-method (ReloadedArraySeed), 10	
altReloadDelayedObjectFunction		
(storeDelayedObject), 16	DelayedMask, 4	
altSaveObject, 17, 18	DelayedMask-class (DelayedMask), 4	
altStoreDelayedObject	DelayedMatrix, 20	
(storeDelayedObject), 16	DelayedSubset, 18	
altStoreDelayedObjectFunction	DelayedUnaryIsoOp, 10, 19	
(storeDelayedObject), 16	denseMatrix, 12	
AmalgamatedArray, 2	dim, DelayedMask-method (DelayedMask), 4	
AmalgamatedArray-class	dimnames, DelayedMask-method	
(AmalgamatedArray), 2	(DelayedMask), 4	
AmalgamatedArraySeed	extract_array,DelayedMask-method	
(AmalgamatedArray), 2	(DelayedMask), 4	
AmalgamatedArraySeed-class	extract_sparse_array,DelayedMask-method	
(AmalgamatedArray), 2	(DelayedMask), 4	
AmalgamatedMatrix-class	extractComponents (AmalgamatedArray), 2	
(AmalgamatedArray), 2	3	
	H5SparseMatrix, 9, 20	
chunkdim, DelayedMask-method	HDF5Array, 9	
(DelayedMask), 4		
	is.na, 5	
cloneDirectory, 11, 12 coerce, AmalgamatedArray, AmalgamatedMatrix-me	is_sparse,DelayedMask-method	
(AmalgamatedArray), 2	(DelayedMask), 4	
coerce, AmalgamatedMatrix, AmalgamatedArray-me	thoddannay	
(AmalgamatedArray), 2	loadArray(readArray), 6	
coerce, ReloadedArray, ReloadedMatrix-method	loadWrapperArray (WrapperArraySeed), 19	
(ReloadedArraySeed), 10	Todowi apper Array (wi apper Array Seed), 19	
coerce, ReloadedMatrix, ReloadedArray-method	matrixClass,AmalgamatedArray-method	
(ReloadedArraySeed), 10	(AmalgamatedArray), 2	
componentNames (AmalgamatedArray), 2	matrixClass,ReloadedArray-method	
createDedupSession, 12	(ReloadedArraySeed), 10	

INDEX 23

path, DelayedMask-method (DelayedMask), 4	stageObject,DelayedMatrix-method
path,ReloadedArraySeed-method	(saveDelayedArray), 13
(ReloadedArraySeed), 10	<pre>stageObject,Matrix-method(saveArray),</pre>
preserveDelayedOperations, 3 , 5	11
	storeDelayedObject, <i>12–14</i> , 16, <i>18</i>
readArray, 6, <i>12</i>	storeDelayedObject,ANY-method
readDelayedArray, 7, 19	(storeDelayedObject), 16
readObject, 6-8, 10	storeDelayedObject,array-method
readSparseMatrix, 8, 15	(storeDelayedObject), 16
recycleHdf5Files, 3, 9	storeDelayedObject,ConstantArraySeed-method
registerReloadDelayedObjectFunction	(storeDelayedObject), 16
(storeDelayedObject), 16	storeDelayedObject,DelayedAbind-method
reloadDelayedObject, 7, 8	(storeDelayedObject), 16
reloadDelayedObject	storeDelayedObject,DelayedAperm-method
(storeDelayedObject), 16	(storeDelayedObject), 16
reloadDelayedObjectFunctionRegistry	storeDelayedObject,DelayedNaryIsoOp-method
(storeDelayedObject), 16	(storeDelayedObject), 16
ReloadedArray, 7, 11, 19	storeDelayedObject,DelayedSetDimnames-method
ReloadedArray (ReloadedArraySeed), 10	(storeDelayedObject), 16
ReloadedArray-class	storeDelayedObject,DelayedSubassign-method
(ReloadedArraySeed), 10	(storeDelayedObject), 16
ReloadedArraySeed, 10, 11	storeDelayedObject,DelayedSubset-method
ReloadedArraySeed-class	(storeDelayedObject), 16
(ReloadedArraySeed), 10	storeDelayedObject,DelayedUnaryIsoOpStack-method
ReloadedMatrix, 8	(storeDelayedObject), 16
ReloadedMatrix-class	storeDelayedObject,DelayedUnaryIsoOpWithArgs-method
(ReloadedArraySeed), 10	(storeDelayedObject), 16
3 //	storeDelayedObject,denseMatrix-method
saveArray, 11	(storeDelayedObject), 16
saveDelayedArray, 13	storeDelayedObject,sparseMatrix-method
saveObject, 6-8, 10, 13	(storeDelayedObject), 16
<pre>saveObject,array-method(saveArray), 11</pre>	storeDelayedObject,SVT_SparseMatrix-method
saveObject,DelayedArray-method	(storeDelayedObject), 16
(saveDelayedArray), 13	(Storebelayedobject), 10
saveObject,denseMatrix-method	type, 5
(saveArray), 11	· 31 · · · · ·
saveObject,ReloadedArray-method	WrapperArray (WrapperArraySeed), 19
(ReloadedArraySeed), 10	WrapperArray-class (WrapperArraySeed),
saveObject,sparseMatrix-method	19
(saveSparseMatrix), 14	WrapperArraySeed, 19
saveObject,SVT_SparseMatrix-method	WrapperArraySeed-class
(saveSparseMatrix), 14	(WrapperArraySeed), 19
saveSparseMatrix, 14	writeSparseMatrix, 20
stageObject, 2, 3, 5, 6, 9	
stageObject, AmalgamatedArray-method	
(AmalgamatedArray), 2	
stageObject, array-method (saveArray), 11	
stageObject, DelayedArray-method	

(saveDelayedArray), 13